

7 Simple Steps to Independence.

Car/Boat/RV Usage

After you are familiar with operating your Inogen One® G4 at home, you are ready to take your Inogen One® G4 on the road.

As always, be a careful driver. DO NOT RUSH OR HURRY to set up your Inogen One® G4 in the car. Take your time and follow the steps outlined to the right. The Inogen One® G4 can operate either by battery or the DC Power Cord while in your car, boat, or RV.

CAUTION

Make sure your Inogen One® G4 is secure before driving.



Plug the DC Power cord into the cigarette lighter or auxiliary DC power supply.



Plug the power output cord into your Inogen One® G4 and position your Inogen One® G4 securely inside the car. Operate your Inogen One® G4 and use normally.



326 Bollay Drive,
Goleta, CA 93117
1-877-4-INOGEN (Toll Free)
1-877-446-6436 (Toll Free)
1-805-562-0515 (International)
info@inogen.net
www.inogen.com

Your Equipment Provider:

Your Physician:

INOGEN ONE^{G4}



INOGEN ONE^{G4}

Now that you have unpacked each box, you have everything you need to start using your Inogen One[®] G4 (IO-400). The set-up of your Inogen One[®] G4 is quick and simple. However, before using the Inogen One[®] G4, please read



Place your Inogen One[®] G4 in a well ventilated location. Ensure the particle filter is in place.



Insert the Inogen One[®] G4 battery by sliding battery into place until the latch returns to the upper position.



Connect the AC input plug to the power supply, connect the power output plug to the Inogen One[®] G4.





Plug the power supply AC cord into the nearest electrical outlet.



Connect the nasal cannula tubing into the metal nozzle on your Inogen One[®] G4. Oxygen will only be delivered while you are breathing through the nasal cannula connected to the



Turn on the Inogen One[®] G4 by pressing the  on/off button.



Set the flow setting to the flow rate prescribed by your clinical professional. The flow setting can be changed by pressing the circular flow setting buttons that are marked with "+" and "-" symbols. Position the nasal cannula on your face and breathe normally through your nose.

