# Inogen Universal Power Supply Instructions for Use

# Inogen Universal Power Supply Instructions for Use

## OVERVIEW

The Inogen Universal Power Supply (BA-107) is used to power the Inogen One concentrator from an AC power source or a DC power source. This allows the user to power the concentrator in the home, in a vehicle or other places where AC or DC power is available.

**NOTE** For use only with the Inogen One Oxygen Concentrator. See the Inogen One Patient or Owner Manual for additional instructions.

## DESCRIPTION


The Inogen One Universal Power Supply is specifically designed for use with the Inogen One Oxygen Concentrator (IO-100). The Universal Power Supply provides the precise current and voltage required to safely power the Inogen One, and is designed to operate from specified AC and DC power sources. When used with AC power sources, the supply automatically adapts to input voltages from 100V to 240V permitting use with most power sources throughout the world.

The Universal Power Supply will charge the Inogen One Battery when used with AC input power or when used with a DC power supply, such as the one found in your car. Due to aircraft power limitations, the Universal Power Supply cannot be used to charge the Inogen One Battery when used on aircraft.

The Universal Power Supply is used with the following components:

- Power Supply with attached power output cable to connect to the Inogen One
- Automobile DC Power Input Cable
- AC Power Input Cable
- Airline DC Power Input Cable

**DC cigarette lighter Power Plug for use in Automobile/RV/Boat/Aircraft**

**Power Output Plug**

**Universal Power Supply**

**AC Input**

**DC Input**

**LED Indicator Light**

**AC Power Plug**

**DC 4-Pin Power Plug for use on an Aircraft**

## **OPERATING INSTRUCTIONS**

Ensure the Power Supply is in a well ventilated location as it relies on air circulation for heat dissipation. It is normal for the Power Supply to be warm during operation.

### **For Use At Home With AC Power**

For operation of the Power Supply using an AC power source, follow these instructions:

1. Connect the AC Input to the Power Supply.
2. Connect the AC Power Plug to the power source and connect the Power Output Plug to the Inogen One. The green LED will be illuminated, indicating the Power Supply has input power.

**NOTE** It is not recommended to use an AC plug adapter with the Universal Power Supply. Use Inogen specified AC Power Input cables for international use.

### **For Use In Automobile/RV/Boat**

For operation of the Power Supply using a DC power source, follow these instructions:

1. Connect the DC Input to the Power Supply.
2. Connect the DC Power Plug (Cigarette Lighter Adapter) to the power source and connect the Power Output Plug to the Inogen One. The green LED will be illuminated, indicating the Power Supply has input power.
3. The plug should insert into the socket without excessive force and securely remain in place. The Cigarette Lighter Adapter comes with a slide switch pre-adjusted to the "narrow" position. This will fit most automobile cigarette lighter sockets. If the Cigarette Lighter Adapter is loose in the socket, adjust the switch to the "wide" position.
4. Ensure that the automobile power socket is adequately fused for the Inogen One power requirements (minimum 15 Amp).

**WARNING** The tip of the Cigarette Adapter Plug becomes HOT when in use. Do not touch the tip immediately after removal from an auto cigarette lighter socket.

**CAUTION** Ensure the automobile power socket is clean of cigarette ash and the adapter plug fits properly, otherwise overheating may occur.

**CAUTION** Do not use the Universal Power Supply with a cigarette plug splitter or with an extension cable. This may cause overheating of the DC power input cable.

**CAUTION** Do not jump start the automobile with the Universal Power Supply connected. This may lead to voltage spikes which could shut down and/or damage the Power Supply.

**CAUTION** When powering the Inogen One in an automobile ensure the engine is running. Operating the device without the engine running may drain the car battery.

## For Use On Airlines

1. Remove the Battery from the Inogen One Oxygen Concentrator. Due to aircraft power limitations, the Universal Power Supply cannot be used to charge the Inogen One Battery when used on aircraft.
2. Connect the appropriate DC Power Plug for the type of available airline power. Check with airline personnel to ensure compatibility.

## General

To remove power, unplug the input cord from its source (ie, AC wall outlet, car cigarette lighter adapter).

**CAUTION** The Universal Power Supply is intended for indoor use only. It is not water resistant.

**CAUTION** Ensure the power supply is powered from only one power source (AC or DC) at any given time. The Universal Power Supply may not operate properly if simultaneously powered using both AC and DC sources.

**CAUTION** Do not disassemble the Universal Power Supply. This may lead to component failure and/or safety risk.

**NOTE** Under certain conditions (see Technical Specifications) the Universal Power Supply may shut down. The green LED will blink or will no longer be illuminated. If this occurs, disconnect the Power Supply for at least 10 seconds and reconnect.

## DC Input Cable Fuse Replacement

The Cigarette Lighter DC Power Plug contains a fuse. If the DC Input Cable is being used with a known good power source and the power supply is not operative (no output power and the green LED is not illuminated), the fuse may need to be replaced.

To replace the fuse, follow these instructions and refer to the photograph below.

1. Remove the tip by unscrewing the retainer. Use a tool if necessary.
2. Remove the retainer, tip and fuse.
3. The spring should remain inside the Cigarette Lighter Adapter housing. If the spring is removed, make sure to replace the spring first before inserting the replacement fuse.
4. Install a replacement fuse, Inogen RP#125 (BUSS MDA-12)\* and reassemble the tip. Ensure the retainer ring is properly seated and tightened.

Cigarette Adapter Plug

Fuse

Retainer Tip

Retainer

**CAUTION** For continued protection against risk of fire, only use fuse specified.

## Cleaning

Wipe the surface using a slightly dampened cloth with mild detergent. Do not immerse the unit. Do not clean the unit with isopropyl alcohol or other solvents.

*\*Contact Inogen Customer Care for assistance.*

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Inogen, Inc.
CELEO Pty Ltd