

# INOGEN ONE<sup>G3</sup> System



Inogen One<sup>®</sup>  
AC Power Supply\*  
Item # BA-301



Inogen One<sup>®</sup> G3  
DC Power Cable\*  
Item # BA-306



Inogen One<sup>®</sup> G3  
Carry Bag\*  
Item # CA-300



Inogen One<sup>®</sup> G3  
Lithium Ion Battery\*  
Up to 4.5 hours run time  
Item #BA-300



Inogen One<sup>®</sup> G3  
Extended Life Lithium Ion Battery\*\*  
Up to 9 hours run time  
Item #BA-316

**Only  
4.8 pounds!**

## The Inogen One<sup>®</sup> G3 System (IS-300):

- Compact and Lightweight
- Noticeably Quiet Operation
- Intelligent Oxygen Delivery
- Designed for 24/7 Use

\*Included with system

\*\*Accessory sold separately

## Oxygen. Anytime. Anywhere<sup>®</sup>

The Inogen One G3 delivers the independence of a portable oxygen concentrator in one of the smallest, lightest, and quietest packages available to the oxygen user today. With flow settings from 1-5, the Inogen One G3 is designed to increase independence for most users 24/7. It can be charged at home or on the go - giving you the freedom of Oxygen. Anytime. Anywhere. <sup>®</sup>

The Inogen One<sup>®</sup> G3 is your one solution, your single solution, for oxygen at home or away.

<b>Weight</b>	4.8 lbs. (includes single battery)
<b>Size</b>	Length: 8.75" Width: 3.0" Height: 8.25" (includes single battery)
<b>Oxygen Flow</b>	Pulse dose delivery system Five flow settings 1, 2, 3, 4, 5
<b>Power</b>	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
<b>Battery</b>	Duration (single battery): Up to 4.5 hours Recharge approximately 2-4 hours with AC or DC power Duration (double battery): Up to 9 hours Recharge approximately 4-8 hours with AC or DC power
<b>Noise Level</b>	39 Decibels*
<b>Warranty</b>	3 year or lifetime limited warranty**
<b>Operation</b>	Simple control functions and easy-to-read LCD display
<b>Use</b>	Designed for 24/7 use

\*At flow setting 2

\*\*See Inogen One<sup>®</sup> Limited Warranty Statement for complete warranty information

Manufactured by:

 inogen

326 Bollay Drive  
Goleta, California 93117  
[www.inogenone.com](http://www.inogenone.com)  
[info@inogen.net](mailto:info@inogen.net)

