

INOGEN ONE^{G3} System



Inogen One[®]
AC Power Supply*
Item # BA-301



Inogen One[®] G3
DC Power Cable*
Item # BA-306



Inogen One[®] G3
Carry Bag*
Item # CA-300



Inogen One[®] G3
Lithium Ion Battery*
Up to 4.7 hours run time
Item #BA-300



Inogen One[®] G3
Extended Life Lithium Ion**
Up to 10 hours run time
Item #BA-316

**Only
4.8 pounds!**

The Inogen One[®] G3 System (IS-300):

- Compact and Lightweight
- Noticeably Quiet Operation
- Intelligent Oxygen Delivery
- Designed for 24/7 Use

*Included with system

**Accessory sold separately

INOGEN ONE^{G3}

Oxygen. Anytime. Anywhere[®]

The Inogen One G3 delivers the independence of a portable oxygen concentrator in one of the smallest, lightest, and quietest packages available to the oxygen user today. With flow settings from 1-5, the Inogen One G3 is designed to increase independence for most users 24/7. It can be charged at home or on the go - giving you the freedom of Oxygen. Anytime. Anywhere. [®]

The Inogen One[®] G3 is your one solution, your single solution, for oxygen at home or away.

Weight	4.8 lbs. (includes single battery)
Size	Length: 8.75" Width: 3.0" Height: 8.25" (includes single battery)
Oxygen Flow	Pulse dose delivery system Five flow settings 1, 2, 3, 4, 5
Power	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
Battery	Duration (single battery): Up to 4.7 hours Recharge approximately 2-4 hours with AC or DC power Duration (double battery): Up to 10 hours Recharge approximately 4-8 hours with AC or DC power
Noise Level	39 Decibels*
Warranty	3 year or lifetime limited warranty**
Operation	Simple control functions and easy-to-read LCD display
Use	Designed for 24/7 use

*At flow setting 2

**See Inogen One[®] Limited Warranty Statement for complete warranty information

Manufactured by:

 inogen

326 Bollay Drive
Goleta, California 93117
www.inogenone.com
info@inogen.net



MKT-0290E DTP