



# INOGEN ONE<sup>®</sup> G4 System



Inogen One<sup>®</sup>  
AC Power Supply\*  
Item # BA-401



Inogen One<sup>®</sup> G4  
DC Power Cable\*  
Item # BA-306



Inogen One<sup>®</sup> G4  
Lithium Ion Battery\*  
Up to 2.7 hours run time  
Item #BA-400



Inogen One<sup>®</sup> G4  
Carry Strap\*  
Item # CA-401



Inogen One<sup>®</sup> G4  
Carry Bag\*  
Item # CA-400



Inogen One<sup>®</sup> G4  
Extended Life Lithium Ion Battery\*\*  
Up to 5 hours run time  
Item #BA-408

## The Inogen One<sup>®</sup> G4 System (IS-400):

- Compact and Lightweight
- Noticeably Quiet Operation
- Intelligent Oxygen Delivery
- Designed for 24/7 Use
- Inogen Connect Bluetooth App

\*Included with system

\*\*Accessory sold separately

## Oxygen. Anytime. Anywhere<sup>®</sup>

The Inogen One G4 delivers the independence of a portable oxygen concentrator in one of the smallest, lightest, and quietest packages available to the oxygen user today. With flow settings from 1-3, the Inogen One G4 is designed to increase independence for most users 24/7. The Inogen One G4 is Bluetooth enabled and compatible with Inogen Connect, the free mobile app. It can be charged at home or on the go - giving you the freedom of Oxygen. Anytime. Anywhere. <sup>®</sup>

The Inogen One<sup>®</sup> G4 is your one solution, your single solution, for oxygen at home or away.

<b>Weight</b>	2.8 lbs. (includes single battery)
<b>Size</b>	Length: 5.91" Width: 2.68" Height: 7.2" (includes single battery)
<b>Oxygen Flow</b>	Pulse dose delivery system Three flow settings 1, 2, 3
<b>Power</b>	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
<b>Battery</b>	Duration (single battery): Up to 2.7 hours Recharge approximately 3 hours with AC or DC power Duration (double battery): Up to 5 hours Recharge approximately 5 hours with AC or DC power
<b>Noise Level</b>	40 dBA*
<b>Warranty</b>	3 year standard warranty**
<b>Operation</b>	Simple control functions and easy-to-read LCD display Inogen Connect mobile app compatibility
<b>Use</b>	Designed for 24/7 use and conforms to all applicable FAA requirements for POC use on board an aircraft

\*At flow setting 2

\*\*See Inogen One<sup>®</sup> Limited Warranty Statement for complete warranty information

Manufactured by:

**inogen**

326 Bollay Drive  
Goleta, California 93117  
[www.inogen.com](http://www.inogen.com)

