



YOUR GUIDE TO RESPIRATORY PORTABILITY



Call **1-800-819-0646** for a free consultation



MOBILITY IS A BREATH AWAY.

If you require and have been prescribed supplemental oxygen therapy, get ready for a welcome breath of fresh air. For over 20 years, Inogen has helped patients worldwide replace bulky oxygen tanks with light weight, portable oxygen concentrator systems that deliver supplemental oxygen with the extra advantage of greater mobility.

From the convenience of our popular Inogen One series to its reliable, pulse-dose based Intelligent Delivery Technology, Inogen helps patients receive the benefits of oxygen therapy while maintaining their active lifestyles. This guide will introduce you to Inogen's line of oxygen therapy solutions. To confirm the right option for you, call **1-800-819-0646** to speak with an oxygen specialist.



Inogen **At Home**



Inogen One **G4**



Inogen One **G5**



BUILDING ON A HISTORY OF SUCCESS.

For over twenty years, Inogen has been more than just a leader in the respiratory care industry - we've become a true ally in the well-being of our customers. Our Inogen One line of portable oxygen concentrators (POCs) have become a trusted solution that over a million patients worldwide depend on to live their lives with peace of mind to leave the home.

Visit [Inogen.com](https://www.inogen.com) for more stories of reclaimed mobility.

“THIS IS THE BEST OXYGEN CONCENTRATOR I COULD HAVE EVER BOUGHT.”

—Wayne W.*

Our team of oxygen specialists is passionate about providing greater mobility to the patients who place their trust in us.

Along with being recognized for medical design and high quality standards, Inogen is an Accredited Homecare Provider. This allows us to extend our service philosophy by offering customers access to expert staff able to understand many facets of respiratory care, including issues relating to loss of mobility, breathlessness and fatigue.



*The opinions expressed herein are independently prepared and for informational purposes only. The views expressed are of the individuals' only. The individual/s featured in this content did not receive any financial compensation for providing their testimonial/s to Inogen, Inc. to use.



Nationwide Accredited
Home care provider



A+ rating



Designed in
the USA.



Highest quality standards
with a Management System
certified to ISO 13485.

SUPPORT TO HELP YOU ENJOY GREATER MOBILITY

Our goal is to ensure easy, affordable access to patients who can benefit from an Inogen portable oxygen system. That begins by working with you to create a customized package based on your lifestyle and needs, and offering a 30-day risk-free trial* to help you discover the benefits of Inogen's portable oxygen concentrators yourself.

Ordering your Inogen system is easy, too. Here are two options to help you regain your mobility in no time:

Purchase

Customers are welcome to directly purchase from us. Just call **1-800-819-0646** to speak with an oxygen specialist who will perform a confidential evaluation, and present appropriate Inogen solutions for you. Benefits include:

- **No ongoing co-pays**
- **A new Inogen system covered by a 3-year warranty****
- **Optional lifetime warranty available for your new Inogen system****

Medicare/Insurance

You may be eligible to have your Inogen system covered by Medicare or insurance.*** Our team will directly work with your doctor, Medicare and insurance provider to identify what benefits may be available to you. Call **1-800-819-0646** for more information.

* Applicable for new purchases only.

** See limited warranty statement on www.inogen.com for complete warranty information.

*** Co-payments and deductibles will apply.



Introducing the Inogen Connect App

Customers who select the Inogen One G4 or G5 can enjoy the extra advantage of our new mobile app. Easily connect your smartphone to your Inogen One to assist with set up, monitor system status, adjust settings and more.*** Ask us for details when you call.

Apple App Store



Google Play App Store



Inogen Connect App



***Inogen Connect is optional and is not required to operate the product. Inogen Connect is available for customers in the United States and Canada.

Visit Inogen.com for more stories of reclaimed mobility.

INOGEN ONE.

Portable oxygen to help you enjoy life.

Inogen One is a portable oxygen concentrator designed for mobility. It replaces heavy oxygen tanks, in a system that can be discretely carried in a shoulder bag. While Inogen One is a pulse dose technology system, it's also light weight. Take it on-the-go with confidence. Charging is also a snap, with the ability to plug into a wall outlet, car DC outlet, or with an external battery charger.

Inogen One systems are available with up to six flow settings to fit your needs. Call us to find the model that is right for you.

G4 carry bag

Mobility Takes Flight

*All Inogen One portable systems
meet FAA requirements for air travel.**

6



Call 1-800-819-0646 to reclaim your mobility.



Inogen One **G5**



Inogen One **G4**

*Based on technical performance test data.

ACCESSORIES TO MAKE YOUR INOGEN ONE PERSONAL.

As an Inogen One customer, you can select from a wide range of accessories. Here are a few of our most popular items.

- AC and DC power cables
- Additional batteries
- External battery chargers
- Backpacks
- Replaceable sieve beds and filters
- Hip pack (G4 only)
- Cart (G5 only)
- Carry bags

Browse our latest range and purchase online at inogen.com

Exact accessory styles and colors may slightly vary by product.

“We always used to travel, and I was afraid those days were over. Thanks to the Inogen One, I’m still able to travel anywhere, even overseas.”
—Elye S.*

Hands-free with the G4 hip pack



*The patient represented in this content is a paid consultant for Inogen, Inc. The opinions and experiences presented here are independently prepared and for informational purposes only.

INOGEN *at* **HOME**

OXYGEN THERAPY TO HELP YOU

Our Inogen At Home system combines 5-liter per minute oxygen delivery with technology that is quiet, letting you watch TV or visit with friends with little noise disturbance.

At just 18 pounds, Inogen At Home is a compact unit you can keep at your bedside, or use it during the day as a companion around the house. All with energy efficiency to help you save over time.

INOGEN *at* **HOME**



- ✓ Only 18 pounds
- ✓ Quiet
- ✓ Low power consumption

8




inogen Call 1-800-819-0646 to reclaim your mobility.

Compact, Quiet & Powerful!



SAVE ON YOUR ELECTRIC BILL EACH YEAR.*

Thanks to its energy technology, Inogen At Home is an energy efficient option for oxygen concentrators used at home. Because it's designed for 24/7 use, the power consumption can significantly add up.




* Calculations are based on published US average hourly kilowatt rates, and assume stationary oxygen concentrator usage of 20 hours per day. A flow setting of 2 liters per minute was used as the Inogen At Home power draw, with the published specifications of a common competitive stationary concentrator used as the standard market concentrator power draw.

**The patient represented in this content is a paid consultant for Inogen, Inc. The opinions and experiences presented here are independently prepared and for informational purposes only.

“It's given me a great deal of freedom. With the Inogen, my quality of life has increased.”

—Claudia H.**

Powerful oxygen therapy in a compact size.



Visit [Inogen.com](https://www.inogen.com) for more stories of reclaimed mobility.

DISCOVER THE RIGHT INOGEN SYSTEM FOR YOU.



Inogen **At Home**



Inogen One **G5**



Inogen One **G4**

WEIGHT	18 lbs	4.7 lbs	2.8 lbs
OXYGEN FLOW (liters/min)	Continuous flow 1-5	Pulse-dose, Intelligent Delivery Technology® 1-6	Pulse-dose, Intelligent Delivery Technology® 1-3
BATTERY LIFE	N/A	Single battery: up to 6.5 hours Double battery: up to 13 hours	Single battery: up to 2.7 hours Double battery: up to 5 hours
POWER SENSING	AC Power Supply 100-240 V, 275 W, 50-60 Hz Auto sensing to allow for worldwide use	AC Power Supply 100-240 V, 240 W, 50-60 Hz Auto sensing to allow for worldwide use	AC Power Supply 100-240 V, 240 W, 50-60 Hz Auto sensing to allow for worldwide use
WARRANTY	3-year or lifetime warranty*	3-year or lifetime warranty*	3-year or lifetime warranty*
USE	Stationary	Stationary, portable, travel	Stationary, portable, travel
SIZE	13" L x 7" W 16.5" High	8.15" L x 3.26" W 7.19" High	5.91" L x 2.68" W 7.2" High
NOISE	40 dBA**	38 dBA**	40 dBA**



OUR CUSTOMERS ALWAYS HAVE THE FINAL WORD.

“ My Inogen system is so much lighter than canisters, and gives me hours of freedom to do what I want. I also love the charging options, and the Inogen Connect app. I only wish I found Inogen sooner. ”

- Jerry L.*

“ As someone with COPD, having Inogen at my side makes me feel safe. Last Christmas, my family took me to Walt Disney World, and I had the best day of my life. Thank you for giving me the ability to live my life again. ”

- Situ G.*

“ I’m beyond satisfied with my Inogen portable and plug-in concentrators. The portable one lets me keep doing the work I enjoy, and walks and light exercise, while the plug-in is so quiet and light weight. It travels with me to visit family, and on trips with my husband. I’m so thankful. ”

- Debra E.*

*The opinions expressed herein are independently prepared and for informational purposes only. The views expressed are of the individuals' only. The individual/s featured in this content did not receive any financial compensation for providing their testimonial/s to Inogen, Inc. to use in our program.

Visit [Inogen.com](https://www.inogen.com) for more stories of reclaimed mobility.

JOIN THE INOGEN FAMILY

- 1 Call 1-800-819-0646 to speak to one of our oxygen therapy specialists.
- 2 Select your preferred Inogen system, and we'll work with your doctor to confirm any necessary prescriptions.
- 3 Enjoy personal set up and use assistance after your Inogen system is delivered.

All major credit cards and checks accepted.



Receive a free Medicare eligibility check.

Find out if your Inogen system is eligible for full or partial coverage under Medicare, or your current insurance.* Simply ask us when you call, and we'll do the work for you, including contact with your doctor as appropriate. You also have the option to rent instead of purchase your system.

Call 1-800-819-0646 for details.

*Co-payments and deductibles may apply.



POWERED TO PROVIDE MEDICAL GRADE OXYGEN ALL DAY EVERYDAY

Inogen portable oxygen therapy has helped nearly 1 million patients in need of supplemental oxygen. It's a responsibility we take with exceptional care, which is why we're so pleased to share this guide with you. We understand that each patient begins their journey by exploring their options and making a decision to trust us as a health partner. When you're ready, our team is here to provide immediate assistance and ensure you find the right Inogen system for your health and lifestyle goals. On behalf of the entire Inogen family, thank you.



“I am very pleased and excited to be part of the Inogen family. My freedom to travel and attend activities has been restored...It is truly light weight and easily managed in the car or at home. Inogen is a blessing to me.”

—Debbie T.*





Follow us on Twitter:
www.twitter.com/InogenInc



Like us on Facebook:
www.facebook.com/inogen



www.instagram.com/InogenInc



www.youtube.com/InogenInc

Discover how people like you have
regained their mobility with Inogen.

www.inogen.com/stories



Call **1-800-819-0646** for a free consultation