



Inogen[®] Portable Oxygen Concentrators meet FAA acceptance criteria for carriage and use on board aircraft.



Inogen BA-516 Battery
Inogen BA-316 Battery

According to Federal Aviation Administration Rules and Regulations (114 CFR Parts 1, 11, 121, et al.), the Inogen[®] BA-516 and BA-316 meet the requirement for the carriage of spare lithium-ion batteries by remaining under 100 Wh, and less than approximately 8 grams.

The Inogen BA-516 Battery and Inogen BA-316 Battery contain two, separate 92.2 Wh battery packs each containing the equivalent of 7.68 grams of lithium-ion.

Adhering to this limitation will facilitate passenger carriage and use of POCs on board aircraft and ensure that there are no restrictions on the number of spare lithium ion batteries less than 100 Wh that can be carried on board the aircraft for full operability of the POC throughout the duration of the flight(s).

For more information, please visit:

<https://www.govinfo.gov/content/pkg/FR-2016-05-24/pdf/2016-11918.pdf>