

# INOGEN ONE<sup>G3</sup> System



Inogen One<sup>®</sup>  
AC Power Supply\*  
Item # BA-301



Inogen One<sup>®</sup> G3  
DC Power Cable\*  
Item # BA-306



Inogen One<sup>®</sup> G3  
Carry Bag\*  
Item # CA-300



Inogen One<sup>®</sup> G3  
Lithium Ion Battery\*  
Up to 4.7 hours run time  
Item #BA-300



Inogen One<sup>®</sup> G3  
Extended Life Lithium Ion\*\*  
Up to 10 hours run time  
Item #BA-316

## The Inogen One<sup>®</sup> G3 System (IS-300):

- Compact and Lightweight
- Noticeably Quiet Operation
- Intelligent Oxygen Delivery
- Designed for 24/7 Use

\*Included with system

\*\*Accessory sold separately

**Only  
4.8 pounds!**

## Oxygen. Anytime. Anywhere<sup>®</sup>

The Inogen One G3 delivers the independence of a portable oxygen concentrator in one of the smallest, lightest, and quietest packages available to the oxygen user today. With flow settings from 1-5, the Inogen One G3 is designed to increase independence for most users 24/7. It can be charged at home or on the go - giving you the freedom of Oxygen. Anytime. Anywhere. <sup>®</sup>

The Inogen One<sup>®</sup> G3 is your one solution, your single solution, for oxygen at home or away.

|                    |   |
|--------------------|---|
| <b>Weight</b>      | 4.8 lbs. (includes single battery)  |
| <b>Size</b>        | Length: 8.75"<br>Width: 3.0"<br>Height: 8.25" (includes single battery)   |
| <b>Oxygen Flow</b> | Pulse dose delivery system<br>Five flow settings 1, 2, 3, 4, 5  |
| <b>Power</b>       | AC Power Supply 100-240V, 50-60Hz<br>(auto sensing to allow worldwide use)<br>DC Power Cable: for mobile use in car   |
| <b>Battery</b>     | Duration (single battery): Up to 4.7 hours<br>Recharge approximately 2-4 hours with AC or DC power<br>Duration (double battery): Up to 10 hours<br>Recharge approximately 4-8 hours with AC or DC power |
| <b>Noise Level</b> | 39 Decibels*  |
| <b>Operation</b>   | Simple control functions and easy-to-read LCD display   |
| <b>Use</b>         | Designed for 24/7 use   |

\*At flow setting 2

\*\*See Inogen One<sup>®</sup> Limited Warranty Statement for complete warranty information

Manufactured by:

**inogen**

326 Bollay Drive  
Goleta, California 93117  
[www.inogen.com](http://www.inogen.com)

