


# INOGEN ONE<sup>G3</sup> System



Inogen One<sup>®</sup>  
AC Power Supply\*  
Item # BA-301



Inogen One<sup>®</sup> G3  
DC Power Cable\*  
Item # BA-306



Inogen One<sup>®</sup> G3  
Carry Bag\*  
Item # CA-300



Inogen One<sup>®</sup> G3  
Lithium Ion Battery\*  
Up to 4.7 hours run time  
Item #BA-300



Inogen One<sup>®</sup> G3  
Extended Life Lithium Ion Battery\*\*  
Up to 10 hours run time  
Item #BA-316

## The Inogen One<sup>®</sup> G3 System (IS-300):

- Highest Capacity at or Below 5 Pounds
- Ultra-Quiet, Reliable Operation
- Clinically Proven for 24/7 Use
- Unique Field-Changeable Sieve Beds

\*Included with system

\*\*Battery options available

Additional accessories available; sold separately

**Only  
4.8 pounds!**

## Oxygen. Anytime. Anywhere<sup>®</sup>

The lightweight Inogen One<sup>®</sup> G3 portable oxygen concentrator has the highest capacity of any portable weighing at or below 5 pounds. Patients appreciate its ultra-quiet and reliable operation for 24/7 use. Sieve beds are easily field changeable with no tools required. Advanced diagnostics take the guesswork out of patient support.

The Inogen One G3 offers the patient one, single solution for oxygen at home or away, day or night.

<b>Weight</b>	4.8 lbs. (includes single battery)
<b>Size</b>	Length: 8.75" Width: 3.0" Height: 8.25" (includes single battery)
<b>Oxygen Flow</b>	Pulse dose delivery system Five flow settings 1, 2, 3, 4, 5
<b>Power</b>	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
<b>Battery</b>	Duration (single battery): Up to 4.7 hours Recharge approximately 2-4 hours with AC or DC power Duration (double battery): Up to 10 hours Recharge approximately 4-8 hours with AC or DC power
<b>Noise Level</b>	39 dBA*
<b>Warranty</b>	3 year standard warranty**
<b>Operation</b>	Simple control functions and easy-to-read LCD display
<b>Use</b>	Designed for 24/7 use

\*At flow setting 2

\*\*See Inogen One<sup>®</sup> Limited Warranty Statement for complete warranty information

Manufactured by:

**inogen**

326 Bollay Drive  
Goleta, California 93117  
[www.inogen.com/hme-providers](http://www.inogen.com/hme-providers)

