

INOGEN *at* **HOME**[®]
Home Oxygen Concentrator



Inogen At Home System (GS-100):

Only 18 Pounds
5 Liters Per Minute, Continuous Flow
Extremely Quiet
Energy Efficient



Innovation in Home Oxygen Therapy

The new Inogen At Home is one of the lightest, quietest, and most energy-efficient 5 liter per minute continuous flow oxygen concentrators on the market today.

At approximately 18 pounds, it is half the weight of some other oxygen concentrators in use today. With a low noise level of 40 decibels at flow setting 2 and very low power consumption levels, the Inogen At Home is an efficient choice for your next home oxygen concentrator.

Weight	18 lbs.
Size	Width: 13.0" Depth: 7.0" Height: 16.5"
Oxygen Flow	Continuous flow system Five flow settings from 1-5 liters per minute
Power Sensing	AC Power Supply 100-240V, 275W, 50-60Hz (auto sensing to allow worldwide use)
Noise Level	40 Decibels at flow setting 2
Warranty	3-year warranty*
Operation	Simple control functions
Use	Designed for stationary use

*See Inogen At Home Limited Warranty Statement for complete warranty information

Manufactured by:

326 Bollay Drive
Goleta, California 93117
www.inogen.com
info@inogen.net

