


INOGEN ONE[®] G4 System



Inogen One[®]
AC Power Supply*
Item # BA-401



Inogen One[®] G4
DC Power Cable*
Item # BA-306



Inogen One[®] G4
Lithium Ion Battery*
Up to 2.7 hours run time
Item #BA-400

**Approximately
half the size of the
Inogen One G3!**



Inogen One[®] G4
Carry Strap*
Item # CA-401

**Only
2.8 pounds!**

The Inogen One[®] G4 System (IS-400):

- Unbelievably Compact and Lightweight
- Ultra-Quiet Operation
- Intelligent Pulse Oxygen Delivery
- Designed for 24 Hour Use

*Included with system

**Accessory sold separately



Inogen One[®] G4
Extended Life Lithium Ion Battery**
Up to 5 hours run time
Item #BA-408



Inogen One[®] G4
Bag*
Item # CA-400

INOGEN ONE^{G4}

Oxygen. Anytime. Anywhere[®]

Once again, Inogen innovation sets a new standard in oxygen delivery.

The Inogen One G4 delivers the independence of a portable oxygen concentrator in one of the smallest, lightest, and quietest packages available to the oxygen user today. With flow settings from 1-3, the Inogen One G4 is designed to dramatically increase independence for most users 24/7. It can be charged at home or on the go - giving you the freedom of Oxygen. Anytime. Anywhere. [®]

The Inogen One[®] G4 is your one solution, your single solution, for oxygen at home or away.

Weight	2.8 lbs. (includes single battery)
Size	Length: 5.91" Width: 2.68" Height: 7.2" (includes single battery)
Oxygen Flow	Pulse dose delivery system Three flow settings 1, 2, 3
Power	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
Battery	Duration (single battery): Up to 2.7 hours Recharge up to 3 hours with AC or DC power Duration (double battery): Up to 5 hours Recharge up to 5 hours with AC or DC power
Noise Level	40 dBA*
Warranty	3 year or lifetime limited warranty**
Operation	Simple control functions and easy-to-read LCD display
Use	Designed for 24/7 use and conforms to all applicable FAA requirements for POC use on board an aircraft

*At flow setting 2

**See Inogen One[®] Limited Warranty Statement for complete warranty information

Manufactured by:



326 Bollay Drive
Goleta, California 93117
www.inogen.com



MKT-0379C