

## **I N O G E N O N E**<sup>G5</sup>

## Oxygen. Anytime. Anywhere®

## Once again, Inogen innovation sets a new standard in oxygen delivery.

The Inogen One G5 offers the most oxygen per pound for a portable oxygen concentrator on the market today. With flow settings from 1-6, the Inogen One G5 is designed to dramatically increase independence for most users 24/7. The Inogen One G5's user interface has been designed with an easy-to-read LCD display and simple control functions. It can be charged at home or on the go - giving you the freedom of Oxygen. Anytime. Anywhere.

Weight 4.7 lbs. (includes single battery) Size Length: 7.19" Width: 3.26" Height: 8.15" (includes single battery) **Oxygen Flow** Pulse dose delivery system Six flow settings 1, 2, 3, 4, 5, 6 **Power** AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car **Battery** Duration (single battery): Up to 6.5 hours Recharge up to 3 hours with AC or DC power Duration (double battery): Up to 13 hours Recharge up to 6 hours with AC or DC power **Noise Level** 38 dBA\* Warranty 3 year or lifetime limited warranty\*\* Operation Simple control functions and easy-to-read LCD display Inogen Connect App Use Designed for 24/7 use and conforms to all applicable FAA requirements for POC use on board an aircraft

The Inogen One G5 is your one solution, your single solution, for oxygen at home or away.

\*At flow setting 2

\*\*See Inogen One® Limited Warranty Statement for complete warranty information

Manufactured by:



