Stationary Oxygen Concentrator

Inogen® Voxi™ 5


Quality Meets Affordability

The Inogen® Voxi™ 5 is an affordable stationary oxygen concentrator solution, with features including:

- Flow Range: 1.0 – 5.0 L/min
- Oxygen Concentration at 5 L/min: 87 – 96%
- Power Supply: 120V, 60Hz
- Input Power: 350VA
- Compact[1] with castors
- Quiet performance: 45 dBA typical[2]
- Audible alarm
- 5-Year expected service life
- 10,000 ft. altitude rating



Catalog Number GS-200-V2



Introducing the

Inogen® Voxi™ 5 Stationary Oxygen Concentrator

The compact and durable Inogen® Voxi™ 5 system offers up to 5-liters per minute continuous flow oxygen delivery from the trusted name in oxygen therapy. Your patients can feel comfort and confidence using the Inogen® Voxi™ 5 at home, whether during the day or at night.

To learn more, call your local Inogen® representative at 855-631-2438 or email us at B2Bsalessupport@inogen.net

SPECIFICATIONS

Catalog Numbers	GS-200-V2 Inogen® Voxi™ 5 Stationary Oxygen Concentrator
Weight	34.2 lbs. (15.5 kg)
Size	Width: 15.35 in. (39 cm) Depth: 9.65 in. (24.5 cm) Height: 19.69 in. (50 cm)
Sound Level	45 dBA[2]
Oxygen Flow	1-5 L/min
Oxygen Concentration	87% - 96% at 1-5 L/min
AC Power Supply	120V - 60Hz, Input Power: 350VA
Expected Service Life	5 years
Accessories	RP-452 Filter RP-453 Firesafe Valve RP-456 Humidifier Bottle (Salter Labs®) RP-457 Filter Cover

1. 2,922 in3/0.048 m3
2. When measured at 2 m (6.6 ft) from device at 3 L/min and 5 L/min

Inogen® Voxi™ 5 Stationary Oxygen Concentrator
INDICATIONS FOR USE: The Inogen® Voxi 5 Stationary Oxygen Concentrator is used on a prescriptive basis by patients requiring supplemental oxygen. It supplies a high concentration of oxygen and is used with a nasal cannula to channel oxygen from the concentrator to the patient. The Inogen® Voxi 5 Stationary Oxygen Concentrator may be used in a professional healthcare facility and home healthcare environment. DO NOT USE IF: This device is not intended to be used in any way other than described in the indications for use. This device is to be used as an oxygen supplement and is not intended to be life sustaining or life supporting. ONLY use this product if the patient is capable of spontaneous breath, able to inhale and exhale without the use of a machine. Not for pediatric use. Not for use by tracheotomized patients. WARNINGS: The device produces enriched oxygen gas, which accelerates combustion. Do not use the device when smoking or near open flame, matches, petroleum, oil, grease, solvents, radiant heaters, aerosols, etc. If you feel ill or uncomfortable, consult your equipment provider and/or your health care provider immediately. If you are unable to communicate discomfort, you may require additional monitoring and or a distributed alarm system to convey the information about the discomfort and or the medical urgency to your responsible caregiver to avoid harm. Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns. To avoid danger of choking or strangulation hazard, keep cords away from children and pets. TALK TO YOUR HEALTH CARE PROVIDER: The oxygen flow setting must be determined and recorded for each patient individually by the prescriber, including the configuration of the device, its parts, and the accessories. It is the responsibility of the patient to periodically reassess the setting(s) of the therapy for effectiveness. The proper placement and positioning of the prongs of the nasal cannula in the nose is critical for oxygen to be delivered.

PM250090 EN_EX_NA | Rx Only. © 2026 Inogen, Inc. | 859 Ward Drive, Suite 200, Goleta, CA 93111
Inogen® and other related product names contained herein are trademarks of Inogen, Inc. The usage of any Inogen, Inc. trademark is strictly forbidden without the prior consent of Inogen, Inc. All other trademarks are trademarks of their respective owners or holders.