inogen

OXYGEN CONCENTRATORS

Your Guide to Portable Oxygen



Call **1-800-819-0646**
for a free consultation

INOGEN®
PORTABLE OXYGEN CONCENTRATORS

Oxygen therapy that moves with you.



Introducing the Inogen Rove Family

If you require and have been prescribed supplemental oxygen therapy, get ready for a welcome breath of fresh air with the Inogen Rove 4™ and Rove 6™.

For over 20 years, Inogen has helped over 1 million patients replace bulky portable oxygen tanks with lightweight, portable oxygen concentrators (POCs) that deliver supplemental oxygen while on the go.

As America's leader in portable oxygen therapy, Inogen helps patients receive the benefits of oxygen therapy while maintaining their active lifestyles. This guide will introduce you to our line of oxygen therapy solutions. **There's only *one* Inogen.**

Call 1-800-819-0646
to speak with an oxygen specialist
and confirm the right option for you.

INOGEN®
PORTABLE OXYGEN CONCENTRATORS

Our compact alternative to tanks.

Portable oxygen concentrators (POCs) provide a safe source of oxygen-enriched air to deliver medical-grade oxygen as prescribed by a physician, freeing patients from managing heavy, cumbersome tanks and planning for refills.

Inogen POCs work by pulling in the air around you, separating the nitrogen, and delivering a concentrated amount of oxygen on each breath.

With reliable, pulse-dose based Intelligent Delivery Technology®, Inogen POCs are powered to provide medical grade oxygen delivery all day, every day. Pulse dose oxygen is delivered based primarily on breathing rate as a sensor determines when a patient's inhale begins to deliver an oxygen pulse at that moment. Patients get the oxygen they need, when they need it.

1 At flow setting 1
2 See limited warranty statement on www.inogen.com for complete warranty information.
3 At flow setting 2
4 When used with a standard battery

Call **1-800-819-0646** to get started



7.50 in.

Inogen Rove 4™

- Only 2.9 lbs. (with standard battery)
- 4 Pulse-dose, Intelligent Delivery Technology® Settings (1-4)
- Up to 3.00 hours run time with standard battery; up to 5.75 hours with extended battery[1]
- 3-year or lifetime warranty[2]
- 39 dBA[3] operating volume
- AC Power Supply 100-240 V, 100 W, 50-60 Hz; Auto sensing to allow for worldwide use
- 5.9" L x 2.7" W



8.14 in.

Inogen Rove 6™

- Only 4.8 lbs (with standard battery)
- 6 Pulse-dose, Intelligent Delivery Technology® Settings (1-6)
- Up to 6.25 hours run time with standard battery; up to 12.75 hours with extended battery[1]
- 3-year or lifetime warranty[2]
- 39 dBA[3] operating volume
- AC Power Supply 100-240 V, 100 W, 50-60 Hz; Auto sensing to allow for worldwide use
- 7.18" L x 3.27" W

INOGEN®
PORTABLE OXYGEN CONCENTRATORS

Accessories for life on-the-go.

Purchase carry bags, extra batteries, power cables, replacement filters, and more accessories to keep you moving at **inogen.com.**



Rove 6 Carry Bag


Rove 4 Extended Batteries



All Inogen portable systems meet FAA acceptance criteria for carriage and use onboard aircraft.*

Call **1-800-819-0646** to get started

* Based on technical performance test data.



** The opinions expressed herein are independently prepared and for informational purposes only. The views expressed are of the individuals' only. The individual featured in this content did not receive any financial compensation for providing their testimonial/s to Inogen, Inc. to use in our program.



“I’m beyond satisfied with my Inogen® portable and plug-in concentrators. The portable one lets me keep doing the work I enjoy, and walks and light exercise, while the plug-in is so quiet and light weight. It travels with me to visit family, and on trips with my husband. I’m so thankful.

Debra E.*

* The opinions expressed herein are independently prepared and for informational purposes only. The views expressed are of the individuals’ only. The individual featured in this content did not receive any financial compensation for providing their testimonial/s to Inogen, Inc. to use in our program.



INOGEN®
STATIONARY OXYGEN CONCENTRATOR

Compact and quiet at-home therapy.

Inogen At Home®

Our Inogen At Home system combines 5-liter per minute continuous flow oxygen delivery with technology that is quiet, letting you watch TV or visit with friends with little noise disturbance.


At just 18 pounds, Inogen At Home is a compact unit you can keep at your bedside, pack in carry on luggage, or use it during the day as a companion around the house. The At Home is energy-efficient to help you save in electric costs over time.

Call 1-800-819-0646
to speak with an oxygen specialist
and confirm the right option for you.

Inogen At Home®

- 18.0 lbs.
- 5 Continuous flow settings (1-5)
- 3-year or lifetime warranty[1]
- 45 dBA[2] operating volume
- AC Power Supply 100-240 V, 275 W, 50-60 Hz; Auto sensing to allow for worldwide use





INOGEN®
STATIONARY OXYGEN CONCENTRATOR

Power & efficiency in a compact size.

Save on your electric bill each year*

Thanks to its energy technology, Inogen At Home® is an energy efficient option for oxygen concentrators used at home. Because it's designed for 24/7 use, the power consumption can significantly add up.

$365.55 per year

Standard Market Stationary Concentrator

$93.73 per year

Inogen At Home Concentrator

* Calculations are based on published US average hourly kilowatt rates, and assume stationary oxygen concentrator usage of 20 hours per day. A flow setting of 2 liters per minute was used as the Inogen At Home® power draw, with the published specifications of a common competitive stationary concentrator used as the standard market concentrator power draw. Actual costs and/or savings may vary.

1 See limited warranty statement on www.inogen.com for complete warranty information.

2 At flow setting 2

Call **1-800-819-0646** to get started

NEXT STEPS

Joining the Inogen Family

- Your doctor provides a prescription saying that you require supplemental oxygen and/or have a severe lung disease.
- Your medical documentation indicates that you are mobile in your home and would benefit from the use of a portable system.
- Skip the tank. Inogen concentrators are available day one of your oxygen therapy when prescribed a POC.
- Call **1-800-819-0646**. We'll check your portable oxygen concentrator eligibility with Medicare, private insurance, or walk you through cash pay options.



inogen®

SCAN THE CODE
to discover how people like you are enjoying their mobility with Inogen®



DESIGNED & ASSEMBLED USA

BBB ACCREDITED BUSINESS A+ rating

bsi. ISO 13485 Medical Devices Quality Management

1-800-819-0646
Call for a free consultation

Inogen, Inc.
859 Ward Drive, Suite 200, Goleta, California 93111
Copyright © 2025 Inogen | EX_EN_USA | PM250036